

# OVERLINE

ITALIA



*Infinity*


Technology Made In Italy 

# TOTAL BODY SCULPT

THE PERFECT SYNERGY FOR A FABULOUS BODY



-  ENDOMASSAGE
-  ELECTROLIPOLYSIS
-  1 MHz ULTRASOUND
-  ELECTROPORATION
-  RADIO FREQUENCY
-  CHROMO FREQUENCY
-  LLLT
-  SHOCK WAVE
-  VACUUM



# Ultra remodelling synergy Cellulite\* Adipose & Atony

A cutting edge appliance that uses a mixture of technologies that are effective in remodelling the silhouette which, thanks to special handpieces, can be used synergetically or individually.

It allows reducing, elasticising, toning, modelling and draining actions thanks to the lipolytic stimulating action of the Low Level Laser (LLLT) which acts directly, in a targeted and non-invasive manner on the cell membrane, favouring the exchange and elimination of fat from the adipose cells.

The appliance is also equipped with Endomassage technology for draining, 1MHz Ultrasound for a slimming action, Electrolipolysis for toning, Radio Frequency for a Collagen stimulating action, Electroporation for a transdermal action, Chromo Frequency for a rebalancing action. Let's not forget that, being a complete technology, it also has special handpieces for the face, to offer a total look solution for any customer.

- ★ REVOLUTIONARY SYNERGY
- ★ LASER REDEFINITION
- ★ MULTIFUNCTION HANDPIECES
- ★ ERGONOMICITY
- ★ REDEFINITION OF THE SILHOUETTE
- ★ MIX OF SPECIFIC RADIO FREQUENCY
- ★ VEHICULATION OF ACTIVE INGREDIENTS
- ★ TONING ACTION
- ★ RELAXATION OPTION
- ★ INTENSIVE FUNCTION
- ★ TRANSDERMAL ACTION

*\*imperfections caused by cellulite.*



## CUSTOMISED PROTOCOLS

Thanks to extensive studies to create this sophisticated technology, a high level software was developed that places the user in the condition to perform any type of treatment on their customer.

It begins with the assessment of the Body Mass Index (BMI), and then proceeds with a test based on simple questions, which helps understand the morphology of the customer being treated. The somatotypes, which have been studied since ancient times, are nothing more than the different body types.

No one of us is identical to anyone else and no one has the absolute characteristics referred to one of the three, very often an individual matches one or more somatotypes. Thanks to this software and its algorithm we are able to establish the dominant one and proceed with the most accurate treatment.

Each somatotype has specific characteristics: for example, imperfections caused by cellulite, imperfections caused by adipose deposits or even the combination of the two. Thanks to this technology we are able to determine the body's specific imperfections and treat them with the suitable protocol. The software is also able to establish the necessary number of sessions required to improve the situation. This was not designed to downplay the professional skills of the Beautician, quite the opposite, we want to support them in achieving reassuring results. If you wish to put yourself to the test without this support, simply choose the free treatment.

***No one has ever offered you more choice than this!***





**SYNERGY ENDO**  
*endomassage - LLLT*  
*chrome frequency - radio frequency*




**BIPOLAR RF**

# TOTAL BO

**TECHNO WAVE**  
*shock wave massage - radio frequency*  
*LLLT - vacuum*






**SYNERGY ULTRA**  
*1MHz ultrasound- electrolipolysis  
chromofrequency - radio frequency  
electroporation*

**MULTIPOLAR RF**

# ODY SCULPT



# TECHNO WAVE

Thanks to its research laboratories, Overline has combined multiple energies into a single movement. The synergetic handpiece, with four specific actions that are developed in unison, is able to achieve **concrete and visible results from the very first session.**



# SYNERGY ENDO

Thanks to studies carried out by the research and development laboratory, Overline has developed multiple treatments in a single action. The synergetic handpiece, with four specific actions that are developed in unison, is able to achieve **concrete and visible results from the very first session.**




# SYNERGY ULTRA

Overline has generated the possibility of performing multiple treatments in a single movement thanks to the technologies featured in the handpiece, multiple treatments in a single movement. The synergetic handpiece, equipped with five specific actions that develop simultaneously, is able to achieve **remarkable results from the very first sessions.**




# RF MANUAL

For a more localised action on particularly hard-to-reach areas, Overline has added two handpieces, with slightly reduced dimensions, to facilitate a more specific and dedicated job




## ELECTROLIPOLYSIS

Electrolipolysis treatment is available among the various options of this technology. The objective of this energy consists in redefining the unwanted imperfections caused by adipose deposits effortlessly and without generating any undesirable effects.



## ULTRASOUND

Ultrasound are waves that travel at a certain speed based on the density of the tissue. They are normally used for common imperfections caused by cellulite and adipose deposits. It is important to remember that the lower the power the greater the penetration.




## ENDOMASSAGE

Specific applicators make it possible to perform rolling, draining and rhythmic massage treatments either simultaneously with other energies or singularly, and this action is indicated for imperfections caused by cellulite.




## VACUUM

The vacuum massage, with other energies, provides remarkable lymphatic draining to clear the body, also improving the connective tissue which can inform the skin structure, making it seem remodelled.




## ELECTROPORATION

Electroporation is a treatment that promotes the penetration of and boosts the vehiculation of dermocosmetics into the deep layers of the skin, allowing them to work directly and immediately.




## SHOCK WAVE

A strong wave of pressure in any elastic medium such as air, water or a solid substance. Thanks to this capacity it is able to act on adipocytes only, remodelling fatty deposits and the imperfections of cellulite.



## RADIO FREQUENCY

RF treatments are used as a method to provide heat to the deepest tissues thanks to the endogenous action of RF, with the aim of increasing collagen production, stimulating circulation and affecting the release of adipose deposits from the tissues.



## LLLT

(Low level laser treatment) Its effects are those that, by generating photo-stimulation, act on toning and firming, reducing oxidative stress and improving cellular metabolism and reducing inflammation.

# REVOLUTIONARY HANDPIECES


# TECHNO WAVE

This extraordinary multifunction handpiece encloses various energies inside of it. It is highly sophisticated handpiece that allows you to achieve never-before-seen results on your customers with a few too many imperfections.

The photos show how it is possible to work on some specific areas. The movements are represented by a grid of circles going upwards, and horizontal during closing, with simple yet effective lines, we will see the body get modelled and adapt under the **TechnoWave** action. The extraordinary synergy allows us to work on areas such as the entire leg, all the way down to the ankle. We can work on the sides of the abdomen to create tension and feel fitter. We can also treat the shoulders for a relaxing action. The objective that we have set with this technology is to resolve the imperfections caused by cellulite and adipose deposits that affect so many people.








## SYNERGY


The underlying principle of this handpiece consists in the combination of multiple energies to achieve a common objective, i.e. to obtain tangible results and surprise customer expectations. Thanks to the commitment of Overline laboratories we have succeeded in inserting ShockWave, Multipolar Radio Frequency, Vacuum and LLLT (Low Level Laser Treatment). The results are very quickly visible thanks to internal\* tests: on a reduced sample of people, we noted a constant reduction in the general measurements and most of all the redefinition of the main imperfections.



### RADIO FREQUENCY MULTIPOLAR


The combination with radiofrequency stimulates the Collagen fibres, making the supporting fibres flexible and extendible.

This makes it possible to achieve results from the very first sessions and within about ten the result will exceed



### LOW LEVEL LASER TREATMENTS

LLLT technology, a laser dedicated to the stimulation of biological processes and to activating the production of the energy that stimulate the muscles.



### SHOCK WAVE + VACUUM


In combination with Vacuum the shock wave acts locally (on the fold created by the vacuum) by performing a specific action of redistribution of the tissues, leaving them smooth and soft. The fold allows us to work carefully and precisely, offering the possibility of working on imperfections of medium-large intensity or small and focalised.

\* Protocols and performance carried out at the Company on a given sample of people.

# REVOLUTIONARY HANDPIECES SYNERGY ENDO



Revolutionary synergy specific for the imperfections of cellulite\*,  
adipose and relaxation of the skin.



Connectival Massage

3

Endomassage

1

Toning Massage  
Endo + RF


2

Supplied with  
3 interchangeable  
accessories

Body treatments will be more efficient and evident thanks to the **revolutionary multifunctional handpiece**. The ergonomics of this handpiece make it easy to grip, without tiring the beautician's hand rather, they can concentrate exclusively on the beauty programme and the planned objectives.

**Synergy means results:** with the precious combination of different energies, the work on the body becomes more targeted and the results will be quickly achieved. Through the action of **Multipolar Radio Frequency** combined with **LLLT** technology that prepare the tissue and stimulate the fibres to actively react, the imperfections of cellulite, adipose and sagging skin will appear considerably improved. **Endomassage** helps complete the treatment, to achieve a silhouette that will be forever changed. While the **Chromo Frequency** re-harmonises the area with the benefit of the vibration of the colours, providing a sensation of well-being and enhancing beauty.

\* imperfections caused by cellulite.



## RADIO FREQUENCY

*Regenerating Action*

Radio Frequency is a non invasive technique that exploits the effect of radio waves to remodel skin tissues, with a particularly effective action on fighting skin sagging on the body and on the face.



## ENDOMASSAGE

*Draining Action*

Specific applicators make it possible to perform rolling, draining and rhythmic massage treatments either simultaneously with other energies or singularly, and this action is indicated for imperfections caused by cellulite.

## LOW LEVEL LASER TREATMENTS

*LLLT Stimulating Action*

Low Level Laser treatment (LLLT) uses a beam similar to a single frequency electromagnetic wave pencil with 980 nm wave length to promote tissue stimulation and an effective body remodelling effect.



## CHROMO FREQUENCY

*Warming/Rebalancing Action*


The Chromo Frequency considers the imperfection to be a disharmonious vibration of the skin cells, which can be rebalanced by exploiting the energy and therefore the wave length of one or more specific colours.



# REVOLUTIONARY HANDPIECES SYNERGY **ULTRA**



Revolutionary synergy specific for the imperfections of cellulite\*,  
adipose and relaxation of the skin.



This revolutionary handpiece is called Ultra because, **thanks to the synergy of various specific actions, it improves the most common imperfections attacking the silhouette of the body, on both women and men.** The combined action of Ultrasound and Electrolipolysis is crucial to remodeling, especially on non-toned areas. The heat of the LLLT technology and the action of the Chromo Frequency stimulate the revitalisation of the fibres. Electroporation is the best tool to help absorb precious active ingredients in the areas where it is most required. A winning synergy of technologies to counter atony and imperfections.

\* imperfections caused by cellulite.

## 1 MHz ULTRASOUND

*Slimming Action*

The Ultrasound treatment is based on the periodic vibrations of particles of a material medium. Vibrations that propagate through the medium, at a finite speed. Ideal in the treatment of imperfections caused by cellulite and adipose.

## ELECTROLIPOLYSIS

*Toning Action*

The objective of Electrolipolysis is to reduce the imperfections of cellulite and fat deposits. Fat is reduced in the selected areas without any effort or discomfort using low voltage electrical current through skin plates. It is a non-invasive treatment. Electrolipolysis exploits the application of micro-currents to effectively treat the imperfections of cellulite\*.

## ELECTROPORATION

*Intensifying Action*

Electroporation is a treatment that promotes the penetration of and boosts the vehiculation of dermocosmetics into the deep layers of the skin, allowing them to work directly and immediately.

## CHROMO FREQUENCY

*Warming/Rebalancing Action*

The Chromo Frequency considers the imperfection to be a disharmonious vibration of the skin cells, which can be rebalanced by exploiting the energy and therefore the vibration of one or more specific colours.

## LOW LEVEL LASER TREATMENTS

*LLLT Stimulating Action*

The LLLT treatment exploits a beam of single frequency electromagnetic waves and predefined wave length to promote tissue stimulation and an effective body remodelling effect.



## HANDPIECES

# RADIO FREQUENCY

to model & define the Body



### MULTIPOLAR RADIO FREQUENCY

*Remodelling action*

Multipolar Radio Frequency is the ideal solution for fighting the signs of ageing. This technology has the benefit of focusing the radius of action in a single point through the 4-pole system. Accordingly it is possible to achieve much high temperatures on the surface and in the deeper layers. The multipolar handpiece also makes it possible to work with lower powers, and therefore slightly longer treatment times, without creating discomfort and/pain for the customer.

### BIPOLAR RADIOFREQUENCY

*Elasticising Action*

The effect of Radio Frequency is based on the transformation of energy into heat, in fact, very often the perceived temperature does not correspond to the temperature in depth precisely because of the endogenous action of this energy. One of the fundamental characteristics of radio frequency is the capacity to give tissues tone thanks to the heat that stimulates the collagen fibres to degenerate and then regenerate stronger than ever. It restores some physiological functions and most of all the treatment is very relaxing and does not involve any post-recovery.



# SMART RADIO FREQUENCY

## SMART HANDPIECES

Smart Radio Frequency Handpiece that helps the beautician perform the Radio Frequency treatment.

## TEMPERATURE READER

Its main feature is that it reads the temperature of the treated area, when the technology is paused, displaying it on the screen of the handpiece and guaranteeing compliance **with the regulations issued by the Health Ministry, with law decree No. 206 which entered into force in January 2016, relative to the rise in temperature during Radio Frequency treatments.**

## LIGHT SIGNAL

Thanks to the evolved software that the equipment was implemented with, the **SMART** handpiece recognises the selected treatment and lights up, identifying which handpiece to use.

## EXECUTION TIME INDICATOR

The time indicator on the display tells the beautician the total energy execution times.

## INTEGRATED LEDS FOR CHROMO WELL-BEING

The integrated LEDs make it possible to perform a Chromo Well-being treatment in synergy, which consists in irradiating areas of the body with the vibration of different coloured lights, to achieve an aesthetic benefit.





## EDUCATIONAL VIDEOS

Within the guided programme, there is also the possibility of viewing some videos that help clarify any doubts and improve dexterity. Cutting edge in terms of service, **Overline Infinity Total Body Sculpt** has always been there alongside professional beauticians to offer excellent results. Standing strong on this quality at **Overline Infinity Total Body Sculpt**, some Educational Videos have been provided to allow the beautician to consult, at any time, the correct operational protocols. Educational videos on the single protocols and manual manoeuvres can be easily consulted on the 10" touchscreen. **Overline** came up with this solution to allow everyone, if they forget minor details, or when new operators are recruited, to use the self-teaching videos in order to achieve the skills required to offer and guarantee the best service.



## PERFECT PERFORMANCE

Protocols and performance carried out at the Company\* in compliance with the suggested work methods, confirm the synergy between the Overline technology and iSol Beauty cosmetics guarantee the maximum efficiency of the treatment.

\* Protocols and performance carried out at the Company on a given sample of people.






## INTUITIVE SOFTWARE

The advanced **Overline Infinity Total Body Sculpt** software can be used for numerous and evolved functions and to control the equipment from the 10" monitor, in a simple, fast and intuitive way. The way it has been conceived, the 10 inch monitor immediately becomes a fundamental tool for implementing the beautician's work. In a few steps, the beautician is able to provide a comprehensive consultation of the customer's specific and physical situation. He/She can examine their customer's images, check programmes and treatments, access data. He/She can also check the situation of what treatment steps to carry out, giving them the possibility of motivating and involving the customer even more. Within the guided programme, it is also possible to view some videos that help clarify any doubts and improve dexterity.





# INFINITELY

# Beauti-

**Overline Infinity** is the innovation choice for the beauty Salon. It is the possibility of infinite customised treatments with the synergy of latest generation technologies. It features an intuitive interface and simple operation. It means effective treatments and Salon customer satisfaction. It involves on-going training support allowing every beautician to make full use of all of the potential offered by the equipment. It provides the safety of using strictly tested, certified systems built in compliance with regulations in force. It has an appealing design of sinuous lines that suits every salon.

**Infinity** is for those who never lose sight of their end goals!

# THE DETAILS

The electronically reclining *Overline Infinity* Monitor makes programming easier and aids control over the work, offering a clear vision of the treatment, also to the customer. A large top drawer and storage compartment for materials and cosmetics, while the shaped shelf on the back of the equipment offers tidy housing for the handpieces.



Storage compartment



Large top drawer



Shaped handpiece shelf

For *Overline Infinity* the concept of customisation is very important. Every technology possesses the advanced *Contactless Card* reading system, where the sequence of each customer's treatment is entered, saving past treatments and those yet to be carried out, so as to optimise the time.



Contactless Card





# 36-MONTH WARRANTY CAREFREE



**Overline Infinity Total Body Sculpt** equipment is guaranteed against materials and labour defects, ex factory, for a period of 36 months. To activate the warranty validity conditions simply fill in every section of the Warranty form inside the **Infinity Total Body Sculpt** package, and send it to **Overline** by and no later than 10 days from the date of purchase.



## INTEGRATED COSMETICS

The duo of *Overline Infinity Total Body Sculpt* and *iSol Beauty* integrated cosmetics is synonymous with concrete solutions for all aesthetic needs and targeted results. *iSol Beauty* cosmetics are designed specifically to work efficiently with the technologies that the equipment features, for the maximum satisfaction of the beautician and an excellent beauty experience for the customer.



- |   |  |  |  |
|---|--|--|--|
| <p>1 Cleansing<br/><b>CLEANSING SOFT MOUSSE</b><br/>ISO.FOAM.150<br/>[150ml]<br/>Dermo-compatible Cleanser<br/>Face&amp;Body</p>  | <p>4 Body treatment<br/><b>RF BODY EMULSION</b><br/>Body Emulsion for RF<br/>ISO.RF.500<br/>[500ml]</p>                    | <p>7 Body treatment<br/><b>RESOLUTION BODY CELL</b><br/>Liquid Crystal Cellulite Serum<br/>for Body Electroporator<br/>ISO.POR.500 [250ml]</p> | <p>10 Body treatment<br/><b>HYDRABODY CELL</b><br/>Cellulite Body Emulsion<br/>ISO.BODY.100<br/>[200ml]</p>    |
| <p>2 Cleansing<br/><b>JOJOBA FACE &amp; BODY SCRUB</b><br/>Gel Scrub With Jojoba Micro-Granules<br/>ISO.SCRUB.200<br/>[250ml]</p> | <p>5 Body treatment<br/><b>ULTRASOUND BODY GEL</b><br/>Innovative gel with Phosphatidylcholine<br/>ISO.GEL.110 [500ml]</p> | <p>8 Body treatment<br/><b>RESOLUTION BODY TONE</b><br/>Liquid Crystal Toning Serum<br/>for Electroporator<br/>ISO.POR.550 [250ml]</p>         | <p>11 Body treatment<br/><b>HYDRABODY FIT</b><br/>Toning Effect Body Emulsion<br/>ISO.BODY.200<br/>[200ml]</p> |
| <p>3</p>  | <p>6 Endomassage &amp; synergy<br/><b>OILGEL SYNERGY</b><br/>ISO.RFOIL.100<br/>[500ml]<br/>RF &amp; Massage</p>            | <p>9 Bo</p>  | <p>12</p>  |

# Always at the service of beauticians

Cutting-edge and Innovation are the ideal terms to define **Overline**. A historical Milanese company, undisputed leader in the production and sale of hi-tech professional cosmetic equipment.

Founded in 1989 and, ever since, it can very simply be described as a “success story”. A company that expresses its business style through ethical commitment and careful and responsible commercial management.

*The literature of the company uses terminology allowed at an international level, as well as being the most commonly used terms to define technologies and features*



# Research & Development

The beating heart of the company, the innovative Research and Development laboratory, is the place where Technologies of high performance and simple use provide effective solutions to the demands of professional beauticians. Various, safer systems that offer multiple functions for face, body and hair removal

## Made in Italy

100% Made in Italy production that offers competitive advantages in terms of performance, safety and service to the customer. Overline technologies are distributed at the finest cosmetic and beauty salons across Italy and in every European country in addition to South America, Mexico, Africa and Asia.

## Certifications

**Overline** is a trademark of Technology S.r.l., a company that constantly pursues the research and development of technologies, to offer the professional beauty market innovative systems that offer visible results, while maintaining the maximum safety of the beautician and customer. Thanks to its great experience, **Technology** possesses important certifications for compliance with regulations and international standards in the production of biomedical, electro-cosmetic and cosmetic equipment. The certifications that make all of the technologies very high level are the following:



**Tüv Rheinland ISO 9001:2015.** This represents the internationally recognised Quality Management System certification that ensures that the devices comply with the management system required by the standard in terms of design, production, sale, development, and research.

**Tüv Rheinland ISO 13485:2016.** The certification related to the Quality Management System of Medical Devices is a standard reference for Quality as part of the regulatory procedures governing the design, production and marketing of Medical Devices.

**Tüv Rheinland ISO 45001:2018.** Relative to the standard for the management of the health and safety of workers.

[www.tuv.com](http://www.tuv.com) ID 9105082302



# OVERLINE

ITALIA



*Infinity*

# ARBIDORA

MEDICAL TRADING CO. L.L.C

عربيدرا لتجارة المعدات الطبية ذ.م.م



CALL +971 50 480 7988

[www.arbidora.com](http://www.arbidora.com)



*Overline reserves the right to make changes to the equipment described, at any time and without prior notice.  
Overline equipment is designed and programmed to obtain results on imperfections without exceeding the allowed cosmetic limit.*



ISO 9001:2015



ISO 13485:2016



ISO 45001:2018

*Overline Infinity* is a Technology S.r.l. trademark